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1. SYSTEM DESCRIPTION 系统描述

Blow-fill-seal (BFS) technology is an automated process by which ampoules are formed, filled, and sealed in a continuous operation. All these procedures are completed in a single time under aseptic conditions, and the whole cycle is completed within 12 to 14 seconds.

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
As the container is formed, filled and sealed on a single machine, the whole process is completed under the protection of aseptic laminar flow air, and all process media (liquid medicine, air, etc.) are sterilized and filtered, so the possibility of the product being polluted by the production environment is very low. During the whole process, the opening time of the container is exposed to the air in the production environment for less than 1.5 seconds, so the BFS aseptic filling process has become the first choice for products that need to be protected from particle and microbial contamination.

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2. PROCESS DESCRIPTION 工艺描述

The following picture show the major phases of BFS machine.

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Process 工艺步骤	Description 描述
<p>Extrusion 挤塑</p> 	<p>High temperature screw extrusion: the thermoplastic material is continuously extruded into tubular shape by high-speed rotating screw, overhanging plastic tube (called parison), and the parison washed by sterile air is through the forming mould. The temperature of plastic hot melt and extrusion is 170~230 degrees, and the pressure is 350bar. It is the aseptic extrusion and heat removal process.</p> <p>高温螺杆挤塑: 用高速旋转螺杆将热塑性材料连续挤压成为管状, 悬垂的塑料管 (称为管坯), 由无菌空气冲洗的管坯穿过成型模具。塑料热融和挤出的温度为 170~230 度, 压力 160-260bar, 是经过实证的无菌挤出和去热源工艺。</p> <p>Parison sealing and cutting: when the parison reaches an appropriate length, the main mould closes and seals the bottom of the parison, the cut-off knife cuts the parison. The top of the parison is positioned by the opening holder, and sterile air continues to protect the opening parison in the mould.</p> <p>管坯封口切断: 当型坯达到适当长度时, 主合模把管坯底部密封, 切刀切断管坯, 管坯顶部被开口保持架定位, 无菌空气继续保护模具中敞开的管坯。</p>
<p>Blowing 吹瓶</p> 	<p>Aseptic air blowing: the mould is quickly transferred to the blow blowing/filling station, and the aseptic air supplied in the grade A air shower protects the blowing, filling and sealing station. Lower the filling needle into the parison until it is sealed against the mould. The plastic material is then molded into containers by vacuum and/or sterile compressed air according to the mould shape.</p> <p>无菌气体吹瓶: 模具被快速转移到吹瓶/灌装工位, A 级风淋装置中供应的无菌空气对吹瓶、灌装和封口工位进行保护。降低灌装针, 使其进入管坯, 直至与模具形成密封。再通过真空和/或无菌压缩空气按照模具的形状, 把塑料原料成型容器。</p>

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Process 工艺步骤	Description 描述
<p>Filling 灌装</p> 	<p>Container quantitative filling: the measured liquid medicine is immediately filled into the container, and sterile air is exhausted. Air and liquid medicine are filtered through a filter before entering molded or forming containers. The pipeline which the liquid medicine passes through should conduct CIP and SIP before filling to ensure the safety of the liquid medicine.</p> <p>容器定量灌装: 通过计量的药液立即被灌装在容器里, 并排出无菌空气。空气和药液在进入已经成型或正在成型的容器之前都经过过滤器过滤。药液经过的管路在灌装之前进行在线 CIP 和在线灭菌 SIP, 确保药液安全性。</p>
<p>Sealing 封口</p> 	<p>Sealing of the ampoule head: after filling, the plastic between the top of the mould and the opening holder is in a state of semi-melting, and then the sealing mould is closed to form the top of the container and make the ampoule sealed.</p> <p>瓶体头部封口: 灌装结束后, 模具顶部和开口保持架之间的塑料处于半融化状态, 接下来密封模具合闭, 形成容器的顶部, 并使瓶子密封。</p>
<p>Outfeeding and cutting 出瓶和切边</p> 	<p>Outfeed, transferring and deflashing: after sealing, the mould is opened, the finished ampoule after forming, filling and sealing is sent out of the machine to start the next cycle.</p> <p>出瓶传送切边: 密封后, 模具打开, 已经完成成型, 灌装并密封的瓶子被送出机器, 接着开始下一个循环。</p> <p>Deflashing system: the deflashing system is placed in the general area. It needs two times of handover from the filling room to general area. The handover area is recommended to be placed in the buffer room to avoid contamination of the clean area caused by cross-domain of different areas.</p> <p>切边系统: 切边系统放置在一般区域, 从灌装间进入一般区经过两次交接完成, 交接区域建议放置在缓冲间, 避免不同区域跨域造成的洁</p>

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The capability of the whole compact line can be adjusted, and the test is subject to the distribution of WFI after the machine running in normal condition.

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
4. APS BLOW-FILL-SEAL MACHINE DESCRIPTION 设备简介


BFS is short for Blow-Fill-Seal machine. It's a fully-automatic machine which can blow the thermoplastic material into container, then finish filling and sealing in continuous operation under A Grade environment.

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4.1. Vacuum infeed 真空上料

Raw material in-feed via the vacuum feeding machine, it can achieve the linkage function with main machine. When the raw material level low, automatically start sucking and infeed, when the raw material level high, automatically stop feeding.

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4.2. Extrusion station 挤出工位

Extruder is a kind of mechanism that converts loose plastic particles into molten plastic in high temperature, high pressure and closed environment by directly driving screw through motor.

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The cast aluminum heater with internal water cooling device and PID control technology are equipped on the screw to provide accurate control of three-stage extrusion temperature for the extruder. The West temperature control module PID controls the temperature of 17 temperature control points, and the temperature control machine precisely controls the plasticizing temperature of the screw segment of the extrusion system.

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4.3. BFS mould shifting and locking 移模锁模部件

4.3.1 Mould carriage movement 模架移动

Adopt the rotating oil cylinder to let the mould carriage part reciprocate the movement in the position of extrusion and filling and sealing through the drive of swing arm.

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Mould carriage movement is controlled by displacement sensor and position control card, which ensures that the repetition accuracy of mould carriage movement is $\pm 0.02\text{mm}$. the mould carriage moves at high speed during the moving process without impact.

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Core control parameters: mould shifting speed 1.2s; parallelism of linear guide of mould shifting is 0.02mm; linear guide relative level 0.02mm, repetition accuracy of mould carriage movement is $\pm 0.02\text{mm}$. , no impact.

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4.3.2 Opening and closing action of master mould 主模开合模动作


The linear hydraulic cylinder is used to drive the mould to close the mould, and the cylinder is controlled by proportional regulating valve. The displacement sensor is used to accurately locate the mould position, and opening and closing action curve is controlled to ensure that the mould moves forward at high speed and closing at low speed.

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Master mould closed 主模闭合

4.3.3 Opening and closing action of sealing mould 封口模开合模动作

The duplex linear hydraulic cylinder is used to drive the mould to close the mould, and the cylinder is controlled by proportional regulating valve. The displacement sensor is used to accurately locate the mould position, and opening and closing action curve is controlled to ensure that the mould moves forward at high speed and closing at low speed.

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The pushing plate action of three-point sliding rail positioning and sealing mould can effectively avoid the phenomenon of raising head for the sealing mould closing caused by the original two-point sliding rail positioning.

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Sealing mould closed 封口模闭合

4.3.3 The feeding action 取料机构

The feeding action adopts servo hydraulic linear oil cylinder. When the feeding rod descends to the middle position, there is the feeding plate. When it descends to the low position, the feeding rod breaks away from the plate. The feeding plate holds the plate for cooling and it is conveyed to the deflashing station.

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4.4. Filling system 灌装系统

The filling system adopts time-pressure method for filling. The servo drives the filling head to move down. The filling system is equipped with buffer tank. The buffer tank adopts Mettler Toledo weighing module to control the liquid level of tank and uses pressure transmitting module to control the filling pressure of tank. It is equipped with high precision and high frequency air replenishing/exhaust valve to ensure the constant pressure of filling buffer.

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◆ **Waste ampoule** identification function: the equipment is equipped with waste ampoule

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identification function, the product is not filled and deflashed due to unqualified forming and will be rejected directly.

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◆**Grade A shower device:** BFS equipment adopts the fan to suck the clean air in the clean room. The air passes through the coarse efficiency filter and HEPA filter and enters the filling area. The air pressure chamber adopts the flow equalizing membrane for supply to reduce the disturbance inside the chamber, which can provide grade A shower for the filling needle and non-sealed ampoule body ensuring the sterility of product. The static pressure chamber of the filling system is equipped with air speed detection and differential pressure detection, at the same time, dust particles, plankton and sedimentation bacteria detection can be selected as required.

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Core control parameters: differential pressure, dust particles, plankton, sedimentation bacteria meet grade A requirements.

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◆Filling head 灌装头

The problem of drug bubbles can be solved better by using a tangential intake pipe. All the control elements are outside and easy to maintain. The diaphragm is controlled by high frequency solenoid valve switching vacuum/compressed air.

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4.5. Mould 模具

Mould: including three parts: ampoule body mould, sealing mould and holder. Each part has separately-adjusted cooling water and vacuum pipeline. The mould is made of imported copper alloy, after heat treatment, it has the advantages of non-deformation, high hardness, high strength, anti-rust and good polishing performance, and long service life. Its own cooling and vacuum channels are connected to the cooling and vacuum system of the equipment to ensure consistent temperature of the products and moulds.

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Mould temperature is controlled by self-excited temperature-controlled water valve, which can automatically adjust the water inlet and outlet of mould to ensure the cooling temperature of product and reduce the generation of waste sheet.

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Core control parameters: mould processing control, flatness 0.015mm, verticality 0.02mm, parallelism 0.02mm, constant temperature control of mould.

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4.8. CIP/SIP system CIP/SIP 系统

BFS equipment has CIP/SIP unit. All valves of CIP/SIP system adopt Gemu pneumatic diaphragm valve and Festo valve terminal automatic control. After the operator has installed the sterilization cup, the system will automatically complete the CIP/SIP program.

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The cold point of the system is designed with temperature monitoring to ensure that all points of the system reach the sterilization temperature. System piping shall be designed in accordance with ASME BPE standards. Horizontal piping shall be designed with slope and drainage piping shall be designed with low point. The piping shall meet 3D requirements.

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System pipeline construction shall be performed according to ISPE GEP[®] standard, and material certification, welding seam report, argon gas protection record, welding sample report and endoscopic inspection report shall be provided.

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4.9. Hydraulic system 液压系统

The hydraulic system includes hydraulic oil tank, overflow valve, accumulator, control valve set, servo-driven low pressure pump head and motor driven high pressure pump head, heat exchanger,

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
etc. The main components are internationally renowned brands. Servo low pressure pump controls the flow and minimizes energy consumption.

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The hydraulic system adopts servo pump, the pump is PID adjusted according to the outlet pressure, and the proportion curve of mould shifting, mould locking, mould opening and mould closing action is controlled, which can effectively reduce the impact of equipment and reduce the energy consumption of the hydraulic pump station by 30%. Nickel-plated oil cylinder is without the risk of paint peeling, more in line with the requirements of pharmaceutical GMP.

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- ◆ Pneumatic system: pneumatic system is divided into compressed air system and aseptic compressed air system, the system consists of filter, valve, pneumatic valve terminal and solenoid valve. Compressed air system is to complete the pneumatic action, sterile air system is to complete the air supply of product forming and CIP/SIP.

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- ◆ Cooling system: the cooling system is divided into two independent loops: chilled water system and cooling water system. The chilled water system is used to cool the reducer, hydraulic system and mould of the extruder system. The cooling water system cools the extruder scrap holder, the extruder parison bottom sealing holder, the take-out holder and the extruder screw barrel.

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- ◆ Vacuum system: the vacuum system includes vacuum pump, vacuum reserve tank, vacuum valve which can be independently controlled by each branch and a vacuum module of analog quantity. Provide vacuum source for product forming.

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All utilities should be tested. When the test value fails to meet the requirements of the set value, the system will alarm and stop.

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4.11. Control system 控制系统

The control system of equipment is mainly composed of the main motion controller, microcomputer touch screen, valve terminal, temperature controller, servo controller, frequency converter, servo card and other control components. The control system design conforms to the requirements of Part 11 and GAMP5.

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